# MC2 & Isolite

Bien-Air<sup>®</sup>  
Dental

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## EMC Precautions

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the user's manual. Essential performance was defined as the maintaining of the visual lightning intensity when powered at 3,5VDC.

### 1. Use of MC2

#### Attention :

The MC2 should not be used adjacent or stacked with other equipment. If adjacent or stacked use is necessary, the MC2 should be observed to verify normal operation in the configuration in which it will be used.

### 2. Use of Accessoires

#### Attention :

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by Bien-Air as replacements parts for internal components, may result in increased emissions or decreased immunity of the MC2.

### 3. Radio Transmitting Equipment

#### Attention :

The MC2 motors comply with the EMC requirements according to IEC 60601-1-2. Radio transmitting equipment, cellular phones, etc. shall not be used in the close proximity of the device since this could influence the performance of the device. Particular precaution must be considered during use of strong emission sources such as High Frequency surgical equipment and similar so that e.g. the HF cables are not routed on or near the device. If in doubt, please contact a qualified technician or Bien-Air Dental.

## Guidance and manufacturer's declaration – electromagnetic emissions

The MC2 motors are intended for use in the electromagnetic environment specified below. The customer or the user of the MC2 motor should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The MC2 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The MC2 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

### Guidance and manufacturer's declaration – electromagnetic immunity

The MC2 motors are intended for use in the electromagnetic environment specified below. The customer or the user of the MC2 motor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±2 kV contact ±4 kV contact ±6 kV contact ±2 kV air ±4 kV air ±8 kV air	±2 kV contact ±4 kV contact ±6 kV contact ±2 kV air ±4 kV air ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Not applicable
Surge IEC 61000-4-5	±0.5 kV line to line ±1 kV line to line ±0.5 kV line to earth ±1 kV line to earth ±2 kV line to earth	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $UT$ (>95 % dip in $UT$ ) for 0,5 cycle 40 % $UT$ (60 % dip in $UT$ ) for 5 cycles 70 % $UT$ (30 % dip in $UT$ ) for 25 cycles <5 % $UT$ (>95 % dip in $UT$ ) for 5 sec	Not applicable	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical industrial location or hospital environment.

NOTE  $UT$  is the a.c. mains voltage prior to application of the test level.

### Guidance and manufacturer's declaration – electromagnetic immunity

The MC2 motors are intended for use in the electromagnetic environment specified below. The customer or the user of the MC2 motors should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	10 V	Portable and mobile RF communications equipment should be used no closer to any part of the MC2, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter  Recommended separation distance  $d = \left[ \frac{3.5}{r_1} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MC2 is used exceeds the applicable RF compliance level above, the motor MC2 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the MC2 motor.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distances between portable and mobile RF communications equipment and the MC2

The MC2 motors are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MC2 motors can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MC2 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of the transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{r_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.035	0.035	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.7
10	1.1	1.1	2.2
100	3.5	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.



NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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This product may be covered by one or more of the following patents:

**EP Europe:** 745358 / 688539 / 948294 / 1145688 / 1563800 / 1563801 / 1675523 / 1753360 **DE Germany:** 29616023.7  
**DK Denmark:** 9600315 **FR France:** 2722972 **CH Switzerland:** 693922 CN China: 100528099 / 100522100 / 100522099 / 100553584 **JP Japan:** 3892485 / 4298933 / 7000419 **US United-States:** 5453008 / 6033220 / 6319003 / 7214060 / 7448870  
**RU Russia:** 2361540 / 2361541 / 2372046

REF 1600681-001 MOT MC2 LED / REF 1600073-001 MOT MC2 IR / REF 1600075-001 MOT MC2 GT AV  
 REF 1600078-001 MOT ISOLITE LK / REF 1600079-001 MOT ISOLITE LI

**Bien-Air Dental SA**  
 Länggasse 60  
 Case postale  
 2500 Bienne 6, Switzerland  
 Tel. +41 (0)32 344 64 64  
 Fax +41 (0)32 344 64 91  
 office@bienair.com

**Bien-Air Italia s.r.l.**  
 Via Vaina 3  
 20122 Milano, Italia  
 Tel. +39 (02) 58 32 12 51/52/54  
 Fax +39 (02) 58 32 12 53  
 ba-i@bienair.com

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**Bien-Air UK Ltd**  
 Arundel House Unit 1 - Ground Floor  
 Amberley Court, Whitworth Road  
 Crawley, West Sussex,  
 RH11 7XL, England  
 Telephone +44 (0)1293 550200  
 Fax: +44 (0)1293 520481  
 ba-uk@bienair.com

**Bien-Air España, SA**  
 Entença, 169 Bajos  
 08029 Barcelona, España  
 Tel. +34 934 25 30 40  
 Fax +34 934 23 98 60  
 ba-e@bienair.com

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### Список патентованных продуктов Bien-Air Dental с защищёнными правами ©:
Technology Service Co. Ltd.