

# ULTRASOUND      ULTRASOUND

## MODULES



**SP 3055**



**SP 4055**

MODE		
Paro	No	0,1 to 1 Watt
Endo	0,6 to 2,3 Watt	0,4 to 2,8 Watt
Scaling	1,7 to 6,6 Watt	2 to 8 Watt
Boost	No	No
FEATURES		
Frequency	29 KHz $\pm$ 0,5 KHz	27 to 32 KHz
Color Coding System	No	No
Output power stage	Single	Single
Signal shape	Sinus	Square
Efficiency	40 %	40 %
TECHNOLOGY		
Feed-Back	No	No
Push-Pull	No	No
Auto tuning	No	Limited
Cruise control	No	No
Light	No	Option

# OUND ULTRASOUND



Xinetic



SP NEUTRON



NEUTRON  
LED

0,1 to 1 Watt	0 to 1 Watt	0 to 1 Watt
0,4 to 2,5 Watt	0,4 to 2,8 Watt	0,4 to 2,8 Watt
1,7 to 8 Watt	2 to 9 Watt	2 to 9 Watt
No	9 to 12 Watt	9 to 12 Watt
27 to 32 KHz	28 to 36 KHz	28 to 36 KHz
Yes	Yes	Yes
Single	Push-Pull	Push-Pull
Square	Sinus	Sinus
40 %	70 %	70 %
No	Yes	Yes
No	Yes	Yes
Limited	Yes	Yes
No	Yes	Yes
Option	Option	Yes